

→ The paramedical stethoscope

CHESTPIECE

Made of anodised aluminium.

Two models:

- ① single head.
- ② dual head.

Color range assorted to tube range:
black, metallic, blue, red, green.

Screw ring

Made of anodized aluminum.

ACOUSTIC DIAPHRAGM

In polyvinyl.

③ ADJUSTABLE BINAURAL


For user's personal comfort
and correct fitting.

④ EXTERNAL SPRING

No sound interference inside the tube.

⑤ «Y» TUBE

Made of PVC
5 colors available:
black, grey, blue, red, green.



The Spengler selection

The Stethocolor is a model especially selected by Spengler. The stethoscope is designed with an external spring that makes for better propagation of sound waves through the tube.

Spengler also proposes a complete range of stethoscopes manufactured in France, from the "cardiology" to the "nursing" model, that covers the full spectrum of auscultation (heart, vascular, lungs abdomen and joints).



Technical specifications

Chestpiece	material	anodized aluminum
	weight	single head: 20 g dual head: 40 g
Pavillon	diameters	large bell (single or dual): Ø 47 mm small bell (dual): Ø 30 mm
Diaphragm	material	polyvinyl
	diameter	large bell (single or dual): Ø 43 mm
Screw ring	material	anodized aluminium
	diameter	large bell (single or dual): Ø 47 mm
Tube		external: Ø 8.5 mm, internal: Ø 4.5 mm
		one piece moulded PVC
		length : 550 mm
Binaurals	external	anodized aluminum adjustable
Spring	external	steel spring
Eartips		polyethylene
Total length		800 mm
Total weight	single head	80 g
	dual head	100 g
Dimensions	zip bag	24 x 14 x 3,5 cm (L x H x D)
References	in zip bag	
Single chestpiece model		
505 810		Stethocolor® black
505 811		Stethocolor® grey
505 812		Stethocolor® blue
505 813		Stethocolor® red
505 816		Stethocolor® green
Dual chestpiece model		
505 860		Stethocolor® black
505 861		Stethocolor® grey
505 862		Stethocolor® blue
505 863		Stethocolor® red
505 866		Stethocolor® green

After-sales service : all spare parts are available from your Spengler dealer.

The Spengler quality

Quality assurance

French manufacturing certified ISO 9001 version 2000, ISO 13485.

CE mark

European directive CEE 93/42 (Class I).
→ The Spengler stethoscopes comply with the stipulations of the guidelines for safety, health, hygiene and protection of the environment.

After sales service

Spengler ensures the maintenance and repairs and trains its registered distributors for these services.

Guarantee

→ According to french law: The Spengler stethoscopes are guaranteed for life against any hidden fault of manufactory.